

PRODUCT DATA SHEET

Clarion™ – 96V Gel Filtration Column Array


Designed Specifically for Use with Automated Liquid Handling Systems

Sorbtech Clarion-96V Column Arrays have been engineered to provide uninterrupted, efficient, automated simultaneous purification of 96 biological samples (proteins, oligonucleotides, or spheroidal nanoparticles) in a fully compliant standard ANSI-SBS format.

Precision packed with SorbaDex-25 or SorbaDex-50 ultrapure dextran gels, the Clarion 96V Column Arrays permit fast removal of small molecules such as buffer salts, dyes, urea, ammonia, biotin, inhibitors, and other small molecular weight impurities.

Clarion-96V Column Arrays Features:

- **A design-specific for uninterrupted, automated well plate purification systems.**
 - **Preloaded with SorbaDex-25 or SorbaDex-50, ultrapure dextran gels.**
 - **Excellent well-to-well, plate-to-plate consistency of performance.**
 - **Process various sample volumes using gravity or light vacuum.**
- **150 – 300µL** – Clarion 96 Gel Filtration Column Array 300-S25M
 - **400µL** – Clarion 96 Gel Filtration Column Array 400-S25M
 - **500µL** – Clarion 96 Gel Filtration Column Array 500-S25M



1. Column Preparation


- a. Carefully remove the desired number of cap strips from the top of the array and then remove the entire bottom sealing foil.
- b. Allow excess column fluid to drain (via gravity) into a suitable waste reservoir. A vacuum of 950 mbar may be used with a manifold to accelerate his process.

2. Column Washing/Equilibration

- a. Wash each column 4 times (approx. 5 mL total) with either deionized water or buffer (use the same buffer for both equilibration and elution)
- b. Allow the wash buffer to drain completely between each aliquot. A vacuum of 950 mbar may be used to speed up the washing process

3. Sample Application

- Load your samples (up to 500 μL) to each column of the array.
- Do not use a vacuum for sample application. If the sample volume is less than 150 μL , add enough wash or equilibration buffer so that the combined volume of each sample equals 150 μL .



4. Elution

- Using the chart below, determine the pre-run and elution volumes specific for your sample size.
- Load the pre-run volume to each column and let it completely enter gel bed. Do not use vacuum.
- Place a collection plate for sample collection under the array.
- Load the correct elution volume to each column and elute the purified sample by gravity.

*Determined using 64 nmol/ml 25-mer oligo in 0.8 M NaCl

Sample Volume	Pre-run Volume	Elution Volume	Oligo Recovery*	Salt Removed
150 μL	200 μL	300 μL	95%	99.9%
200 μL	150 μL	350 μL	94%	99.4%
250 μL	100 μL	400 μL	96%	99.1%
300 μL	0 μL	500 μL	95%	96.2%